

Blue BHD Self-Sealing Autoclave Bags



Typical Application & Characteristics

Recommended for Steam Sterilisation of Pharmaceutical Production Equipment

- Designed for Steam Sterilization
- Incremental improvement compared to cutting and wrapping with BioShield™ and autoclave tape
- Avoid Sterilization Re-Validation
- Steam Sterilization Indicator
- Eliminates need to cut material and the use of particulate generating autoclave tape
- Fast and Consistent
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Resists tearing and punctures
- 2-Ply for Good Microbial Barrier
- Manufactured and packaged in ISO 6 and 8 cleanrooms
- Certificates of Conformance sent with shipment

Material Specifications

PROPERTY	UNITS	BHD WRAPPING MATERIAL
Basis Weight	gm/m ²	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	litres/100sq. cm	110
Suter Water Head (Pressure Test)	mm	400
Dry Mullen Burst	grams/sq. cm	2,150
Flammability	seconds	8

Sizes Available

Comprehensive range of over 20 sizes from 6" x 12" to 46" x 48".
All size requests, within reason, can be catered for.
Samples available on request.

Call 01270 660429 www.batherton-services.co.uk

Email Addresses: enquiries@batherton-services.co.uk Company Number: 8610746 VAT Number: GB 712 892 429

Registered Address: Batherton Services, Stapeley House, London Road, CW5 7JW



ISO 9001
Certificate No. GB2004565